icumedical
human connections

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24 August, 2010

Linda Gioia,
Nurse Manager
Haematology Oncology Australia - Gold Coast
Level 9, Queen Street
QLD, 4215

Re: Administration of Platelets using Plum A+ Infusion Device.

Dear Ms Gioia,

Thank you for your inquiry regarding the administration of platelets using the Plum® Series (devices including the Plum A+™, Plum® 1.6, Plum XL® Series).

Two *in vivo* studies were performed on blood delivery using the Plum® infusion systems. Dr. Toby Simon^{2,3} performed a red cell study at the University of New Mexico. In the final report, Dr. Simon noted the following:

"No detectable damage to stored red blood cells occurred using the equipment tested."

Adding 60 to 100 ml of 0.9% Sodium Chloride Injection, USP to 300 ml of packed red cells can facilitate flow rates to approximate that of whole blood.

Dr. Ed Snyder¹ performed evaluation of platelet integrity at Yale-New Haven Hospital. The authors concluded:

"...these infusion systems **do not significantly injure or activate platelets**", and that they are "...efficacious for transfusing platelet concentrates to thrombocytopenic patients."

I hope that you find this information helpful. Please contact the Medical Information department on 1300 046 774, if you wish to discuss this matter further.

References:

1. Snyder EL, Rinder HM, and Napychank P: In vitro and in vivo evaluation of platelet transfusions administered through an electromechanical infusion pump. *Am J Clin Pathol* 94(1):77-80, (July) 1990.
2. Simon T: Personal communication and report.

Hospira Pty Ltd
ABN 13 107 058 328
Level 3, 390 St Kilda Rd
Melbourne Victoria
Australia 3004
T: 1300 046 774
F: +61 3 9868 0111
www.hospira.com.au

Hospira Pty Ltd is bound by the National Privacy Principles of the Privacy Act 1988 and other applicable privacy laws. We have to provide our Medical Information services or to otherwise respond to your queries. We may disclose your information to our companies or contractors for the same purposes, or if laws require disclosure. For example, laws may require us to provide our information to authorities such as the Therapeutic Goods Administration. Hospira strives to ensure the safe and appropriate use of all products going forward to the Quality Use of Medicines, it may be necessary as a consequence of your enquiry for an appropriate following your receipt of any advice from Hospira, in order to close the matter or provide further assistance. You have the right to access information under the Privacy Act and other laws. You can do this by phoning us on 1300 046 774. You can also find out our practices at www.hospira.com.au or by contacting us.

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